

# Chem-Impex International, Inc.

---

**Tel:** (630) 766-2112

**E-mail:** sales@chemimpex.com

**Shipping and Correspondence:**

935 Dillon Drive

Wood Dale, IL 60191

**Fax:** (630) 766-2218

**Web site:** www.chemimpex.com


**Manufacturing site:**

825 Dillon Drive

Wood Dale, IL 60191

---

## Product Specifications



|                          |                                                                                               |
|--------------------------|-----------------------------------------------------------------------------------------------|
| <b>Catalogue Number</b>  | 34193                                                                                         |
| <b>Product</b>           | <b>Ethyl phenylacetate</b><br>Phenylacetic acid ethyl ester<br>Benzeneacetic acid ethyl ester |
| <b>CAS Number</b>        | 101-97-3                                                                                      |
| <b>Molecular Formula</b> | C <sub>10</sub> H <sub>12</sub> O <sub>2</sub>                                                |
| <b>Molecular Weight</b>  | 164.20                                                                                        |

---

|                         |                                      |
|-------------------------|--------------------------------------|
| <b>Appearance</b>       | Colorless to nearly colorless liquid |
| <b>Odor</b>             | Sweet, honey odor                    |
| <b>Flash Point</b>      | > 93 °C                              |
| <b>Specific Gravity</b> | 1.027 - 1.032 at 25 °C               |
| <b>Refractive Index</b> | n <sub>20/D</sub> 1.496 - 1.500      |
| <b>Anion</b>            | Halogens : To pass test              |
| <b>Acid value</b>       | 0 - 1 mg KOH/g                       |
| <b>Purity by GC</b>     | 98 - 100%                            |
| <b>Grade</b>            | Natural, FG, KOSHER, FEMA 2452       |
| <b>Storage</b>          | Store at RT                          |

**For Laboratory Use Only**